

EndoNaut®

- PAD -

Navigation, safety
and optimization.



- ✓ Less X-ray and contrast doses
- ✓ Optimized stent size
- ✓ Secured device deployment with augmented roadmaps



EndoNaut®

Designed to upgrade your OR at a lower cost

1 Easy generation of angiography panorama

- Clear and comprehensive view of lower limb






No OR / IT / HR constraints

Just plug EndoNaut® to the C-arm video output.
It's ready for use immediately!




2 Delimitation of lesions

- Lesion cues and length measurement



3 Real-time panorama-based roadmapping

- Lowers the need of contrast agent and X-rays
- Saves time during operations
- Compensates for table and patient movements
- Optimizes balloon / stent length sizes for better outcomes

 www.Therenva.com
 contact@therenva.com
 +33 972 52 29 20
 4 rue Jaurès - 35000 Rennes - FRANCE



EndoNaut® obtained FDA 510(k) Substantial Equivalence determination and is a CE-marked (class IIb) medical device designed in France by Therenva SAS. Not reimbursed. EndoNaut® is a registered trademark of Therenva SAS. Before use, please read carefully instructions contained in the manual or label of the device provided to healthcare professionals. Document only for healthcare professionals. Update 2019/09/02. Non contractual document - ©2019 Therenva SAS. All rights reserved.