



Full range of portable intubation endoscopes for complete airway management



Outstanding performance characteristics

The PENTAX Medical portable intubation endoscopes are cordless and ultra-compact, for quick and reliable intubation techniques and complete airway management. This system concept is ideal for use in emergencies or the outpatient department, satisfying the requirements of anaestheologists in a variety of clinical situations.

Maximum working comfort

- Up to 160° tip deflection, a tight bending radius andwide angle of view, simplifying orientation and the positioning of endotracheal tubes.
- Ergonomically designed, lightweight control body enhances user performance.
- PENTAX Medical LSV (Light Source Versatility) Solutions.
 Three light source options: Cordless battery powered
 LSV module; Self-contained can also be A/C adapter light source module; All models can be connected directly to a conventional light source.
- PENTAX Medical CSV (Closed Suction Valve), a mechanism for accurate aspiration capability.

Superb image quality

- PENTAX Medical Superfine Technology delivers bright and clear, high resolution images for superior observations.
- Upgradeable to video imaging through connection to the PENTAX Medical USB Camera Kit.

PENTAX Medical quality & service programme

- Delivers cost effective products of exceptional quality and durability.
- Provides industry leading service plans, consultancy, financing, technical services and training.

Portable Fibre Intubation Endoscopes

Туре	Ultra- Slim	Slim	Standard	Standard	Large Channel
Model	FI-7 RBS	FI-9 RBS	FI-10 RBS	FI-13 RBS	FI-16 RBS
Field of view [°]	95 90		95		
Depth of field [mm]	3–50				
Tip deflection up/down [°]	130/130			160/130	
Insertion tube [ø mm]	2.4	3.1	3.5	4.2	5.2
Distal end [ø mm]	2.4	3.0	3.4	4.1	5.1
Min. instrument channel [ø mm]	-	1.15	1.31	1.72	2.55
Working length [mm]			600		

In the interest of technical progress, specifications may change without notice.



Julius-Vosseler-Straße 104 . 22527 Hamburg . Germany Tel.: +49 40 / 5 61 92 - 0 . Fax: +49 40 / 5 60 42 13 E-mail: medical@pentax.de . www.pentaxmedical.com

