

Beskrivelse	Insertion Tube ø (mm)	Working Channel ø (mm)	Working Length (mm)	View Angle (°)
Bronkoskoper				
EB15-J10 Video Bronchoscope J10 Series	5,2	2	600	120
EB19-J10 Video Bronchoscope J10 Series	6,4	2,8	600	120
FB-18V Fiber Bronchoscope V Series	6	2,7	600	120
FB-19TV Fiber Bronchoscope V Series	6,2	3	600	120
FB-18RBS Fiber Bedside Bronchoscope	6	2,55	600	100
FB-15RBS Fiber Bedside Bronchoscope	4,9	1,95	600	100
FB-15V Fiber Bronchoscope V Series	4,9	2,1	600	120
FB-8V Fiber Bronchoscope V Series	2,8	1,15	600	100