

Beskrivelse	Insertion Tube ø (mm)	Working Channel ø (mm)	Working Length (mm)	View Angle (°)
Naso-Pharyngo-Laryngoskoper (ØNH)				
VNL9-CP VIVIDEO Video Naso-Pharyngo- Laryngoscope CP Series	3,3	0	300	90
VNL8-J10 HD Video J10 Series Naso-Pharyngo-Laryngoscope	2,9	0	300	80
VNL11-J10 HD Video J10 Series Naso-Pharyngo-Laryngoscope	3,6	0	300	80
VNL15-J10 HD Video J10 Series Naso-Pharyngo-Laryngoscope	4,9	2	300	80