

you won't believe what you see



 **EASY MOVING**

The most powerful x-ray mobile system in the market
now with telescopic column for an increased visibility

 **SEDECAL**

YOUR OEM PARTNER



a world of benefits



EASY MOVING

Fully counterbalanced telescopic column improve effortless positioning without draining power of the batteries.

EASY MOVING

right before your eyes

The new mobile Easy Moving with telescopic column provides all the technologies attributes used for SEDECAL mobiles with new innovations thinking in the patient and the professional: power, reliability, high quality imaging in digital X-ray, coverage and now also extreme visibility.

Thanks to the telescopic column the mobile can access everywhere easily and permits the professional watch anything around, while the x-ray unit is moving.



From excelent to perfect

THE MOST POWERFUL IN THE INDUSTRY

Max.
50 kW


Easy Moving is the most powerful high frequency generator mobile in the market, with a range from 20 Kw to 50 Kw Nominal power. The proven reliability of these generators assures **constant performing during the whole day**.

MAXIMUM MOBILITY AND COVERAGE



Sedecal Easy Moving has been designed to perform in any specific location and patient phase: patient on bed, first aid and emergency, orthopedics, pediatric and operating room. **Thanks to the technologies "Easy driven, easy rotation and easy extension", the level of mobility and coverage of patient is maximum.**

DESIGNED FOR AN ADVANCED VISIBILITY



The new generation of Sedecal mobiles with telescopic column solves any kind of obstacle that users and professional could find while they are driving or parking into the clinic. The new unit provides **maximised vision thanks to the telescopic column and the ergonomic design**. Nurses and doctors will be able to access any location at the hospital in a comfortable and safety way.

HIGHLY CONFIGURABLE FOR TOP FLEXIBILITY



This mobile is ready to **combine with any flat panel detector (wireless or detachable)** and also includes fully integrated digital acquisition software.

Provided with 17" Touch Screen display to match the highest industry standards

 **SEDECAL**

YOUR OEM PARTNER

THE MOST POWERFUL IN THE INDUSTRY



MAXIMUM MOBILITY AND COVERAGE



HIGHLY CONFIGURABLE FOR TOP FLEXIBILITY



DESIGNED FOR AN ADVANCED VISIBILITY



SEDECAL

Calle Pelaya 13 - Pol. Ind. Río de Janeiro
28110 Algete (Madrid - SPAIN)
Central Offices:
Tel: (+34) 91 628 05 44
Fax: (+34) 902 190 385
Email: info@sedecal.com

SEDECAL USA Inc

230 Lexington Drive
Buffalo Grove, Illinois 60089
USA
Tel: 1 (847) 394 - 6960
Fax: 1 (847) 394 - 6966
Email: sales@sedecalusa.com



YOUR OEM PARTNER