


When you need a benchtop, self-contained cabinet X-ray to provide sharp, detailed images.



## The XPERT<sup>®</sup> 20 Benchtop Cabinet X-ray System

The XPERT 20 is a compact, user-friendly X-ray imaging solution. It can be utilized for small sample applications including seed counting and analysis, non-destructive testing, and small animal imaging (i.e., zebrafish and single mouse). The system's benchtop configuration and versatile nature provides you with reliable high-resolution images in seconds from the comfort of your lab.




X-ray image of cardiac valves

### XPERT<sup>®</sup> 20 System Features:

- X-ray imaging system designed for placement directly in the lab.
- Up to 1.3x magnification.
- Seed analyzer software available.
- Easy to operate, intuitive system.
- Desktop or mobile platform option for easy transportation.

# The XPERT<sup>®</sup> 20

Benchtop Cabinet X-ray System



X-ray image of bee larvae




X-ray image of a Zebrafish

## Specifications

Detector Size	3" x 5" (7 cm x 12 cm)
Energy Range	20-50 kV
Tube Current	1.0 mA
Size (W x D x H)	11" x 12" x 20" (28 x 30.5 x 51 cm)
Weight	70 lbs. (32 kg)
Mobile Cart	Optional

**KUBTEC**<sup>®</sup>  
S C I E N T I F I C

M1260B-0521 (11/22) © KUB Technologies Inc. 2022 All rights reserved. KUBTEC, XPERT and associated logos are trademarks and/or registered trademarks of KUB Technologies, Inc., and/or its subsidiaries in the United States and/or other countries. Made in USA 

Scan for more  
product options

